

# Simply Tranquil®



ReadyWear  
Sound Generator  
**ITC**

## DEVICE FEATURES

- Small, In-The-Ear Design
- #10 Battery
- Soft Silicone Body
- Tranquil Circuitry
- Hypo Allergenic
- Volume Control
- 0 dB (HL) VC Floor

## CONSUMER APPEAL

- Cosmetic Appeal
- Economical Price
- No Ear Impression
- Same Day Fitting


## Soothing sounds of the Tranquil all day and all night long.

Simply Tranquil is the preferred choice of tinnitus patients who need a sound generator they can use all day and night!. The smallest, and most comfortable low-level sound generator on the market, it offers natural sound open canal format in an ultra small design. It's soft all-soft body provides a comfortable and stable fit. The open venting permits the application of a stable, broad-band noise without occlusion. Due to the soft, flexible nature of the body, the device seats better in the ear and stays in place throughout the night. The silicone body is hypo allergenic and constantly flexes with jaw movement. Simply tranquil feaures the the popular Tranquil™ brand sound generator.

## DATA SPECIFICATIONS

Source.....	0.0 dB
Peak .....	52.5 dB
Peak Frequency .....	5,900 Hz
RMS OUT.....	63.3 dBSPL
Noise Reduction .....	16 dB
Battery Drain .....	0.39 mA
Battery Life (10A).....	205 Hrs.

## FREQUENCY RESPONSE



**General Hearing**  
Instruments, Inc.

P.O. Box 23748  
New Orleans, LA 70183-0748  
800.824-3021 Fax (504) 733-3767